

CS 3700



Safety, Regulatory and Technical Specifications User Guide

Notice

The Safety, Regulatory, and Technical Specifications User Guide includes information on the safety instructions, regulatory information, and technical specifications of the devices. We recommend that you thoroughly familiarize yourself with this guide to make the most effective use of your system.

© Carestream Dental LLC, 2022. The information in this document is subject to change. Neither Carestream Dental LLC nor any of its subsidiaries shall be liable for errors contained herein or for incidental damages in conjunction with the furnishing, performance, or use of this material. No part of this publication may be reproduced without the permission of Carestream Dental, LLC.

This document is originally written in English.

All trademarks and registered trademarks are the property of their respective holders.

The CS 3700 is intended for professional use only.

U.S. Federal law restricts this device to sale by or on the order of a dentist.

If any serious incident occurs in relation to the device, the user must report it to Carestream Dental and to the competent authority of its Member State in the European Union.

Manual Name: *CS 3700 Safety, Regulatory, and Technical Specifications User Guide*

Part Number: TA1307

Revision Number: 03

Print Date: 2022-02

The CS 3700 complies with Medical Device Regulation (EU) 2017/745 and Medical Devices Regulations 2002 (SI618) as subsequently amended by the EU Exit Regulations of 2019 (SI 791) and 2020 (SI 1478).



Contents

Chapter 1	Indications for Use	1
Safety Information	Clinical Benefits and Performance	
	Characteristics	1
	Conventions in This Guide	2
	Warnings and Safety	
	Instructions	3
	Cleaning, Disinfecting, Sterilizing	6
	Cleaning and Disinfecting the Scanner	6
	Cleaning and Sterilizing the Scanner Tips	8
	Precautions Before Use	12
	Cleaning, Disinfecting, and Sterilizing	12
	Visually Inspecting the Scanner for Damage	12
	Visually Inspecting the Scanner Tips for Damage	12
	Marking and Labeling Symbols	13
	Label Locations	14
	CS 3700 Labels	14
Chapter 2	General Regulatory Information	17
Regulatory	Guidance and Manufacturer's	
Information	Declarations	20
	Compliance with International Regulations	24
Chapter 3	Factory	25
Technical	Manufacturer	25
Specifications	Model.	25
	CS 3700 Technical Specifications	26
	Length of Cables Supplied with the Unit	27
	CS 3700 Environmental Requirements	27
	Computer System Requirements	28

Chapter 4	Manufacturer's Address	29
Contact Information	Authorized Representatives	29
	List of Importers for European Union According to the MDR 2017/745.	30

1 Safety Information

Indications for Use

The CS 3700 is a digital optical scanning device used to record the topographic characteristics of teeth or dental impressions in three dimensions. The resulting topographic impressions are intended for use in the computer-aided design and manufacturing of dental restorative prosthetic devices, dental implant prosthetic devices, and orthodontic models.

The CS 3700 scanner is intended to be used by a dental healthcare professional, for pediatric and adult patients requiring a dental examination.

Clinical Benefits and Performance Characteristics

Carestream Dental's intraoral scanners benefit a dental practice by enabling practitioners to acquire digital impressions with the quality and accuracy required for digital CAD/CAM dental applications. The actual performance of the device is dependent on the user's training and operating execution. The user is solely responsible for the accuracy, completeness, and adequacy of the acquired data.

Conventions in This Guide

The following special messages emphasize information or indicate potential risks to personnel or equipment.



WARNING: Warns you to avoid injury to yourself or others by following the safety instructions precisely.



Caution: Alerts you to a condition that might cause serious damage.



Important: Alerts you to a condition that might cause problems.



Note: Emphasizes important information.



Tip: Provides extra information and hints.

Warnings and Safety Instructions



DANGER OF ELECTRIC SHOCK

This is an electrical unit. Do NOT expose it to water spray. Such action can cause an electric shock or a malfunction of the unit.



IMPORTANT: All known residual risks, contraindications, or undesirable side effects are listed in this guide. If any serious incident occurs in relation to the device, you must report it to Carestream Dental and to the competent authority of your Member State in the European Union.



WARNINGS

CS 3700:

- You **MUST** read and understand this safety information before using the scanner.
- This scanner shall only be used inside hospitals and other professional healthcare facilities and **MUST NOT** be used near high frequency surgical equipment and the RF shielded room of an ME System for magnetic resonance imaging, where the intensity of electromagnetic disturbance is high.
- Before using the scanner, check the outer surfaces of the unit and any accessories to ensure there are no rough surfaces, sharp edges, or protrusions which may cause a safety hazard.
- You are responsible for the operation and maintenance of the scanner. You **MUST** have training to use the scanner.
- **DO NOT** place objects within the field of operation of the unit.
- When the unit is not in use, ensure that the scanner is turned **OFF**.
- **DO NOT** use the scanner in conjunction with oxygen-rich environments. This unit is not intended for use with flammable anesthetics or flammable agents.
- **DO NOT** pull or twist the cable.

- DO NOT drop the scanner or the accessories.
- DO NOT heat sterilize the scanner.
- DO NOT expose the scanner to water spray or submerge it in water or disinfectant.
- DO NOT expose the scanner to high vibrations.
- DO NOT expose the scanner to ultraviolet radiation directly. The scanner is not designed for ultraviolet disinfection.
- DO NOT stare at the LED emission window.
- When the tip is removed, DO NOT touch the heater.
- When the tip is removed, install the rubber sleeve to protect the scanner lens window.
- DO NOT remove the cover of any scanner components. The scanner contains no user-serviceable parts. For any repairs, contact a qualified Carestream Dental service technician.
- DO NOT replace the cables provided with the scanner with other cables. Doing so may damage the scanner and adversely affect the safety protection and EMC performance of the scanner.
- DO NOT replace the power adapter provided with the scanner with any other power adapter. Substitutes may not provide the required protection against electric shocks and other safety hazards.
- Any other equipment not complying with IEC 60601 shall be kept at least 1.83 meters away from the patient.
- If the equipment is faulty, turn it OFF, display an “Out of Service” notice, and contact a qualified Carestream Dental service technician.
- Using components, accessories, cables and spare parts other than those specified or provided by the manufacturer of this equipment may impair the safety protection of the scanner and may result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- No modification of this equipment is allowed.
- Additional multiple outlet strips or extension cords should not be connected to the system.
- The maximum temperature of the applied part may reach to 43 °C; to avoid overheating, do not use it for extended periods.
- To power off the device, push the power button for 3 seconds. To isolate the device from mains supply on all poles, unplug the adapter from the mains power outlet.
- DO NOT maintain or service this equipment while it is in use with the patient.

- Always position the unit in such a way that it is easy to disconnect the adapter from the mains power outlet.
- Connection of the PEMS (Programmable Electrical Medical System) to an IT-NETWORK that includes other equipment could result in risks to patients, operators, or third parties. The responsible organization should identify, analyze, evaluate, and control these risks.

Computer:

- DO NOT place the computer and the peripheral equipment connected to it in the immediate vicinity of the patient. Leave at least 1.83 meters distance between the patient and the equipment.
- The scanner is only intended to be connected to a computer that has agency approval according to the latest edition of applicable safety and EMC standards. Connecting the scanner to other equipment may be hazardous.
- See the installation guide for your computer for information about the data processing system, computer, and screen. Leave a sufficient amount of clear space around the computer to ensure that it is properly ventilated.
- Position the screen to avoid light reflections from internal or external lighting for maximum image quality and visual comfort.

Disposal:



This equipment contains certain materials and chemical compounds incidental to the manufacture of electrical and electronic equipment, and improper "end-of-life" disposal of such equipment can result in environmental contamination. Therefore, this equipment should not be disposed of as ordinary household waste but should instead be delivered to a designated electrical and electronic waste disposal or recycling center. For further information on disposing of electrical and electronic waste, contact the cognizant authority within the local jurisdiction.

Dispose of the scanner tips according to standard operating procedures or local regulations for the disposal of contaminated medical waste. For additional scanner tips, contact your dealer.

Cleaning, Disinfecting, Sterilizing

Cleaning and Disinfecting the Scanner



WARNINGS

- Read and follow the warnings and personal protection instructions provided in the Safety Data Sheet (SDS) for the disinfectant used to process the scanner for reuse.
- You must wear gloves while cleaning and disinfecting the scanner.
- The scanner must be disinfected with an U.S. Environmental Protection Agency (EPA) registered or CE marked intermediate-level disinfectant solution with tuberculocidal activity between patients.
- **DO NOT** use a disinfectant containing phenolics or iodophors; doing so will damage the surface coating of the scanner.
- **Never** put the scanner in an autoclave device or immerse it in water or the disinfectant solution.
- **Excessive fluids can damage the scanner. Do not use cotton, cloth, or tissues soaked with disinfectant to disinfect the scanner.**

Cleaning the Scanner

If the scanner is visibly contaminated with blood and/or body fluids, you must clean it before disinfecting it.

To clean the scanner, follow these steps:

- 1 Dampen (**do not soak**) a lint-free cloth with lukewarm water.
- 2 Remove the blood and/or body fluids with the dampened lint-free cloth.

Disinfecting the Scanner

After each patient, the scanner must be thoroughly disinfected.

To adequately disinfect the scanner, follow the disinfectant manufacturer's instructions for the appropriate contact time.



Important: If the scanner is visibly soiled, it must be thoroughly cleaned prior to disinfecting. See [“Cleaning the Scanner.”](#)

To disinfect the scanner, follow these steps:

- 1 Remove the reusable tip.
- 2 Remove all visible soil (see [“Cleaning the Scanner”](#)).
- 3 Use a commercially prepared intermediate level disinfectant wipe. Follow the manufacturer's instructions for contact time.
Approved disinfectant wipes: Mikrozid AF Jumbo Wipes, CaviWipes, Oxivir Tb Wipes, Clorox Healthcare Bleach Germicidal Wipes, PDI Sani-Cloth Bleach Germicidal Wipes.



WARNING: Using a disinfectant that has not been approved may cause damage to the scanner.

- 4 Thoroughly wipe all surfaces of the scanner.



WARNING: Do not rinse.

- 5 Allow to air dry.
- 6 After the scanner has dried, use a clean, lint-free cloth dampened with water to remove residual disinfectant from the surface of the scanner.

Cleaning and Sterilizing the Scanner Tips

Scanner tips received from the manufacturer are NOT sterilized. You must sterilize the tips before the first use.

The removable scanner tips are autoclavable up to 20 cycles. After 20 cycles, discard the tip. **If you limit the exposure time at 134°C to not more than 4 minutes, you can autoclave the tip up to 60 cycles.**



WARNINGS

- **Wear gloves when handling a contaminated scanner tip.**
- **Read and follow the warnings and personal protection instructions provided in the manufacturer's SDS for the detergent used to clean the scanner tip prior to sterilization.**
- **Do not soak the scanner tips in disinfectant overnight.**
- **Dry the scanner tips thoroughly before mounting onto the scanner.**
- **Do not use an ultrasonic cleaning machine to clean the scanner tips.**

Manually Cleaning the Scanner Tips

To manually clean the scanner tips, follow these steps:

- 1 Rinse excess soil from the tip.
- 2 Using a soft brush, apply an enzymatic detergent solution (e.g., Metrex EmPower) to all surfaces.
- 3 Rinse under clean, running water.
- 4 Inspect the tip. If the tip is not clean, repeat the steps.
- 5 Use a lens tissue or lint-free cloth to remove any dust from the mirror in the tip.

Cleaning the Scanner Tips in an Automatic Washer or Disinfector

To clean the scanner tips in an automatic washer or disinfector, follow these steps:

- 1 Rinse excess soil from the tip.
- 2 Using a soft brush, apply an enzymatic detergent solution (e.g., Metrex EmPower) to all surfaces.
- 3 Load the tip into the washer/disinfector equipment.
- 4 Run the cycle per the equipment manufacturer's instructions.
- 5 If the machine does not have an automatic rinse cycle, rinse thoroughly to remove detergent residues by immersing in clean water.
- 6 Use a lens tissue or lint-free cloth to remove any dust from the mirror in the tip.


Sterilizing the Scanner Tips



Note: The removable scanner tips can be autoclaved for up to 20 cycles. After 20 cycles, discard the tip.
If you limit the exposure time at 134°C to not more than 4 minutes, you can autoclave the tip up to 60 cycles.

To sterilize the cleaned scanner tips, follow these steps:

- 1 Fold a 2x2-inch (5x5cm) 4-ply sterile non-woven gauze sponge in half.



- 2 Using a cotton swab, carefully insert the folded gauze into the tip window. Ensure that the mirror is fully covered by the gauze. To facilitate removal, leave an edge of gauze outside the tip window.



- 3 Place the tip in an FDA-cleared or CE-marked sealed sterilization pouch. The pouch should be sealed air-tight. Use either a self-adhesive pouch or a heat-sealed pouch.

- 4 Place the tips in a steam autoclave for the following times:

Pre-Vacuum Autoclave		
Exposure Time at 132°C	Exposure Time at 134°C	Minimum Drying Time
Minimum 4 Minutes	Minimum 3 Minutes	20-30 Minutes

Gravity Autoclave		
Exposure Time at 132°C	Exposure Time at 134°C	Minimum Drying Time
Minimum 15 Minutes	Minimum 10 Minutes	15-30 Minutes



Important: Do NOT exceed 134°C.



Important: Do not exceed 18 minutes of exposure time.



Important: Never autoclave a tip that does not have gauze covering the mirror and is not wrapped, as this will leave stains on the mirror that cannot be removed.

Precautions Before Use

Perform the following activities on your scanner and accessories before use.

Cleaning, Disinfecting, and Sterilizing

To ensure maximum hygienic safety for the patient and to minimize the risk of cross-contamination, carefully perform the following maintenance activities on your scanner and accessories. After each patient:

- Clean and disinfect the scanner. See [“Cleaning and Disinfecting the Scanner” on page 6.](#)
- Clean and sterilize the scanner tip. See [“Cleaning and Sterilizing the Scanner Tips” on page 8.](#)

In the event that you see poor scan quality or an unclear video preview in the software, clean the tip mirror and the scanner's lens window using a microfiber cleaning swab, applying ethanol that is free of impurities.

Visually Inspecting the Scanner for Damage

Visually inspect the scanner for damage or signs of deterioration by doing the following:

- Inspect the scanner's lens window.
- Inspect around the scanner buttons and the cable.

If damage is noted, do not use the scanner and contact your representative.

Visually Inspecting the Scanner Tips for Damage

Visually inspect the scanner tips for signs of deterioration by doing the following:

- Verify that the tip is not damaged and its components are not detached.
- Verify that the tip mirror does not have any smudges or scratches on it. This can affect image quality.












If deterioration is noted, replace the tip.



WARNINGS

- **The lens window on the scanner is a delicate optical component. Mount the rubber sleeve to protect the lens window from damage and dirt when the scanner is not in use.**
- **The mirror in the tip is a delicate optical component. Its clean and undamaged surface is critical to scan quality.**

Marking and Labeling Symbols

	Type BF applied part symbol classification in accordance with IEC 60601 standards.
	Class II equipment
	<p>In the European Union, this symbol indicates: DO NOT discard this product in a trash receptacle; use an appropriate recovery and recycling facility.</p> <p>Contact your local sales representative for additional information on the collection and recovery programs available for this product.</p>
	Manufacturer's address
	Manufactured date
	<p>CAUTION:</p> <p>Consult accompanying documentation.</p>
	Refer to instruction manual/booklet.
	<p>Caution: Hot Surface</p> <p>Surface can be hot and should not be touched.</p>
	Direct current
	Medical device
	Name of the European authorized representative and address of the registered place of business.

Label Locations

CS 3700 Labels

The following figures illustrate the label locations of the CS 3700.

Figure 1 CS 3700 Box Label




Figure 2 CS 3700 Scanner Label





Figure 3 CS 3700 Heater Label



Figure 4 CS 3700 Power Supply Label



2 Regulatory Information

General Regulatory Information

Compliance with European and International Standards	
EN 60601-1 / IEC 60601-1	Medical Electrical Equipment, Part 1: General requirements for basic safety and essential performance
EN 60601-1-2 / IEC 60601-1-2	Medical Electrical Equipment, Part 1-2: General requirements for basic safety and essential performance – Collateral Standard: Electromagnetic disturbances – Requirements and tests
IEC 60601-2-18	Medical Electrical Equipment, Part 2-18: Particular requirements for the basic safety and essential performance of endoscopic equipment
EN 62471 / IEC 62471	Photobiological safety of lamps and lamp systems: Equipment classification, requirements, and User's Guide
EN ISO 17664	Processing of health care products – Information to be provided by the medical device manufacturer for the processing of medical devices
EN 60601-1-6 / IEC 60601-1-6	Medical Electrical Equipment, Part 1-6: General requirements for basic safety and essential performance - Collateral Standard: Usability
EN/IEC 62366-1	Medical Devices – Part 1: Application of usability engineering to medical devices
EN 62304 / IEC 62304	Medical device software - Software life cycle processes

Compliance with European and International Standards	
EN ISO 10993-1	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process
EN ISO 14971	Medical devices - Application of risk management to medical devices
EN ISO 15223-1	Medical devices - Symbols to be used with medical device labels, labeling and information to be supplied - Part 1: General requirements
EN 1041	Information supplied by the manufacturer of medical devices
CAN/CSA-C22.2 No. 60601-1	Medical Electrical Equipment - Part 1: General requirements for basic safety and essential performance
ANSI/AAMI ES60601-1	Medical Electrical Equipment - Part 1: General requirements for basic safety and essential performance

Classification in Accordance with EN/IEC 60601-1	
Type of protection against electric shock	Class II equipment
Degree of protection against electric shock	Type BF Applied Part
Degree of protection against harmful ingress of water	IPX0 Note: After the tip is installed correctly, the tip portion of the device is IPX1.
Mode of operation	Continuous operation
Flammable anesthetics	Not suitable for use in the presence of flammable anesthetics or a mixture of flammable anesthetics with air or oxygen or nitrous oxide.

Conformity with EN/IEC 60601-1-2

IEC 60601-1-2: 2014 EMC requirements and tests, Medical Electrical Equipment including CISPR 11:2009+A1:2010 Group 1, Class B.



Electromagnetic Compatibility Precautions

Medical electrical equipment requires special precautions regarding electromagnetic compatibility (EMC). Medical equipment must be installed and put into service according to the EMC information provided in this documentation.

Other equipment can interfere with communications with the CS 3700, even if the equipment complies with CISPR emissions requirements.

Warning: Portable and mobile RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the CS 3700, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Guidance and Manufacturer's Declarations

**Guidance and Manufacturer's Declaration -
Electromagnetic Immunity (IEC 60601-1-2)**

The CS 3700 is intended for use in the electromagnetic environment specified below. The customer or user of the CS 3700 should assure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment - Guidance
RF Emissions CISPR 11	Group 1	The CS 3700 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class B	The CS 3700 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic Emissions IEC 61000-3-2	Class A	
Voltage Fluctuations/ Flicker Emissions IEC 61000-3-3	Complies	

Electromagnetic Immunity for Equipment and Systems
Fully Compliant with IEC 60601-1-2: 2014

The CS 3700 is intended for use in the electromagnetic environment specified below. The customer or the user of the CS 3700 should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines	±2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV line to line	±1 kV line to line	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions, and voltage variations on power supply input lines IEC 61000-4-11	0% U_T ; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315° 0% U_T ; 1 cycle and 70% U_T ; 25/30 ^a cycles Single phase: at 0° 0% U_T ; 250/300 ^a cycles	0% U_T ; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315° 0% U_T ; 1 cycle and 70% U_T ; 25/30 ^a cycles Single phase: at 0° 0% U_T ; 250/300 ^a cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the CS 3700 requires continued operation during power mains interruptions, it is recommended that the CS 3700 be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

a) e.g., 10/12 means 10 periods at 50 Hz or 12 periods at 60 Hz

NOTE: U_T is the a.c. mains voltage prior to application of the test level.

Guidance and Manufacturer's Declaration - Electromagnetic Immunity (IEC 60601-1-2)

The CS 3700 is intended for use in the electromagnetic environment specified below. The customer or the user of the CS 3700 should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Electromagnetic Environment - Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz 6 Vrms in ISM bands between 150 kHz and 80 MHz ^a	Environment of a professional healthcare facility.
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz	WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the CS 3700 including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

NOTE: Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the CS 3700 is used exceeds the applicable RF compliance level above, the CS 3700 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the CS 3700.

a The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz.

Guidance and Manufacturer's Declaration - Electromagnetic Immunity

For the immunity to proximity fields from RF wireless communications equipment, the CS 3700 is compliant with the test levels specified below, according to IEC60601-1-2 standard. The customer or user of the CS 3700 should assure that it is used in such an environment.

Test Frequency (MHz)	Band (MHz)	Immunity Test Levels
385	380-390	Pulse modulation 18Hz, 27V/m
450	430-470	FM, ± 5 kHz deviation, 1 kHz sine, 28V/m
710		
745	704-787	Pulse modulation 217Hz, 9V/m
780		
810		
870	800-960	Pulse modulation 18Hz, 28V/m
930		
1720		
1845	1700-1990	Pulse modulation 217Hz, 28V/m
1970		
2450	2400-2570	Pulse modulation 217Hz, 28V/m
5240		
5500	5100-5800	Pulse modulation 217Hz, 9V/m
5785		

Compliance with International Regulations

The CS 3700 complies with the following regulations:

- Medical Device Regulation (EU) 2017/745
- FDA Center for Devices & Radiological Health CDRH - Title 21 CFR 872.3661 (USA)
- Medical Devices Regulations (Canada)
- Directive 2011/65/EU on the Restriction Of the use of certain Hazardous Substances in electrical and electronic equipment (ROHS), as amended by Directive (EU) 2015/863
- Medical Devices Regulations 2002 (SI618) as subsequently amended by the EU Exit Regulations of 2019 (SI 791) and 2020 (SI 1478).

3

Technical Specifications

Factory

Rayco (Shanghai) Medical Products Company Limited
Building 7, No. 1510 Chuanqiao Road
China (Shanghai) Pilot Free Trade Zone
201206 Shanghai
PEOPLE'S REPUBLIC OF CHINA

Manufacturer







Carestream Dental LLC
3625 Cumberland Boulevard, Suite 700
Atlanta, GA USA 30339

Model




CS 3700

CS 3700 Technical Specifications

Components	Technical Specifications
Sensor technology	CMOS
Illumination	LED: Amber, Blue, Green
Field of view	13 x 13 mm 13 x 7 mm (posterior tip)
Depth of field	-2 to +12 mm
Anti-fogging technology	Actively heated tip, guaranteed non-fogging operation when used intraorally
Cable length	2.7 m (1.9 m + 0.8 m)
Digital connection	USB 2.0 High Speed
Dimensions without cable	218 x 36 x 58 mm (with normal/side tip)

Components	Technical Specifications
Weight	316 g (excluding scanner cable and power box)
Handpiece	Input: 12 V  2A
Power Box	75 x 21 x 21mm Input: 12 V  2A Output: 12 V  2A
Adapter	Model UES24LCP-120200SPA: Input: 100-240V ~ 50/60Hz, 500mA Output: 12.0V  2.0A

Length of Cables Supplied with the Unit

Illustration of Part	Part Name	Length of Cable (m)
	Scanner	1.9 m
	Power Box	0.8 m
	AC Adapter	1.8 m

CS 3700 Environmental Requirements

Components	Environmental Requirements
Operating Temperature	+5 ~ 30 °C
Transportation and Storage Temperature	-10 ~ 60 °C
Operating Relative Humidity	10 – 85% RH
Transportation and Storage Relative Humidity	10 – 95% RH
Operating Atmospheric Pressure	700 – 1060 hPa
Transportation and Storage Atmospheric Pressure	600 – 1060 hPa

Computer System Requirements

If necessary, you must update your computer system configuration.

Item	Recommended	Minimum
CPU	Laptop: Intel Core i7-7700HQ, Quad CPU, 2.8 GHz Desktop: Intel Core i7-7700K, Quad CPU, 4.2 GHz	Laptop: Intel Core i7-4700QM, Quad CPU, 2.4 GHz Desktop: Intel Core i7-3770, Quad CPU, 3.4 GHz
RAM	16 GB RAM	16 GB RAM
Monitor	Standard CRT/LCD monitor with screen resolution of 1920 X 1080	Standard CRT/LCD monitor with screen resolution of 1440 X 900
Operating system	Windows 10 Professional (64 bit)	Windows 10 Professional (64 bit)
USB port	USB 2.0 high speed port	USB 2.0 high speed port
Video card	Laptop: NVIDIA GeForce GTX 1050 Ti or Quadro P3000 or similar Desktop: NVIDIA GeForce GTX 1050 Ti or similar	Laptop: NVIDIA GeForce GTX 860M or Quadro K3100M or similar Desktop: NVIDIA GeForce GTX 760 or similar
Video card driver	Support OpenGL 4.3 and OpenCL 1.1	Support OpenGL 4.3 and OpenCL 1.1

The computer and its screen should be situated in or close to the operating area, in the visual field of the practitioner when using the CS 3700.



Important: It is **MANDATORY** to check that your system configuration is compatible with the computer system requirements for the CS 3700 software.



Note: Always use Microsoft Windows Update to ensure that the latest security patches are correctly installed.

4

Contact Information

Manufacturer's Address



Carestream Dental LLC

3625 Cumberland Boulevard, Suite 700
Atlanta, GA USA 30339

Authorized Representatives

Authorized Representative in the European Community



Trophy
4, Rue F. Pelloutier
Croissy-Beaubourg
77435 Marne-la-Vallée, Cedex 2
France

UK Responsible Person

Carestream Dental Ltd.
Samantha Bush
Wiltron House, Rutherford Close
Stevenage
Hertfordshire, SG1 2EF
United Kingdom

Authorized Representative in Brazil

Carestream Dental Brasil Eireli
Rua Romualdo Davoli, 65
1º Andar, Sala 01 - São José dos Campos
São Paulo - Brazil
CEP (Zip code): 12238-577


List of Importers for European Union According to the MDR 2017/745

Carestream Dental France SAS
4 Rue F. Pelloutier, Croissy-Beaubourg
77435 Marne-la-Vallée Cedex 2, France

Carestream Dental Germany GmbH
Hedelfinger Str. 60
70327 Stuttgart, Germany

Carestream Dental Spain, S.L.U.
Paseo de la Castellana, 79
Madrid 28046, España

Carestream Dental Italy S.r.l.
Via Mario Idiojmi 3/3
Assago 20090 (MI), Italia

The background of the page features a large blue area with a subtle gradient and a yellow diagonal stripe running from the bottom left towards the top right. A solid yellow horizontal bar is positioned at the bottom of the page.

For more information, visit: carestreamdental.com
To give documentation feedback, visit:
carestreamdental.com/documentationfeedback